

## qPCR Cycler Programs for Venor® GeM qOneStep Kit

The following protocols were created based on internal experience and customer reports. Minerva Biolabs does not warrant or assume responsibility for the performance of these protocols.

### Programming the LightCycler® 2.0

#### Program 1: Pre-incubation

Cycles	1
Analysis Mode	None
<b>Temperature Targets</b>	<b>Segment 1</b>
Target Temperature [°C]	95
Incubation time [min]	2:00
Temperature Transition Rate [°C/s]	20.0
Secondary Target Temperature [°C]	0
Step Size [°C]	0.0
Step Delay [Cycles]	0
Acquisition Mode	None

#### Program 2: Amplification

Cycles	45		
Analysis Mode	Quantification		
<b>Temperature Targets</b>	<b>Segment 1</b>	<b>Segment 2</b>	<b>Segment 3</b>
Target Temperature [°C]	95	55	60
Incubation time [s]	15	30	45
Temperature Transition Rate [°C/s]	20.0	20.0	20.0
Secondary Target Temperature [°C]	0	0	0
Step Size [°C]	0.0	0.0	0.0
Step Delay [Cycles]	0	0	0
Acquisition Mode	None		
Single	None		

#### Program 3: Cooling

Cycles	1
Analysis Mode	None
<b>Temperature Targets</b>	<b>Segment 1</b>
Target Temperature [°C]	40
Incubation time [s]	30
Temperature Transition Rate [°C/s]	20.0
Secondary Target Temperature [°C]	0
Step Size [°C]	0.0
Step Delay [Cycles]	0
Acquisition Mode	None

#### Result Reading:

- Select the fluorescence channels 1 and 2
- Click on *Quantification* to generate amplification plots and Ct-values
- The threshold will be generated automatically.
- Samples showing no significant increase in the amplification plot can be considered as negative.

## qPCR Cycler Programs for Venor® GeM qOneStep Kit

---

### Programming of ABI Prism® 7500

#### Detector Settings:

Target Probe:	Reporter - FAM™ / Quencher - none
Internal Control Probe:	Reporter - HEX™ / Quencher - none

The “ROX Reference” function needs to be disabled, as no ROX™ dye is included in the mix. Activate both detectors for each well.

Measure fluorescence during annealing.

#### Program Step 1: Pre-incubation

Setting	Hold
Temperature	95 °C
Incubation time	3:00 min

#### Program Step 2: Amplification

Cycles	45
Setting	Cycle
Denaturing	95 °C for 30 sec
Annealing	55 °C for 30 sec & data reading
Extension	60 °C for 45 sec

#### Result Reading:


- Enter the following basic settings at the right task bar:
  - Data: Delta RN vs. Cycle
  - Detector: FAM™ and HEX™
  - Line Colour: Well Colour
- Open a new window with for the *Graph settings* by clicking the right mouse button
  - Select the following setting and confirm with ok:
  - Real Time Settings: Linear
  - Y-Axis Post Run Settings: Linear and Auto Scale
  - X-Axis Post Run Settings: Auto Scale
  - Display Options: 2
- Initiate the calculation of the Ct-values and the graph generation by clicking on Analyze within the report window.
- Pull the threshold line into the graph. Adapt the threshold line to the initial linear section of the positive control reaction.
- Samples showing no Ct-value can be considered as negative.

## qPCR Cycler Programs for Venor® GeM qOneStep Kit


### Programming the ABI StepOne™ / StepOne Plus™

Instrument setup using a template file:

- 1 If your device is using the StepOne™ Software v2.3, download the template file “Venor\_template\_StepOnePlus.edt” and/or “Venor\_template\_StepOne.edt” from our web page ([www.minerva-biolabs.com/en/detection/venor-gem-qep](http://www.minerva-biolabs.com/en/detection/venor-gem-qep)) and save it to a folder on your hard drive or your network.
- 2 Start the StepOne™ Software




- 3 For setting up a new experiment using setup information from a template click on the Template button. A new window will pop up.




**qPCR Cycler Programs for Venor® GeM qOneStep Kit**


- Navigate to the folder containing the template and open the appropriate template file (Venor\_template\_StepOnePlus.edt for StepOnePlus 96 well instrument and Venor\_template\_StepOne.edt for StepOne 48 well instrument, respectively.)



- Presented parameters in the sections Experiment Properties, Plate Setup and Run Method will be loaded.




- In the Plate Setup section further sample names can be defined. (Negative control and Positive control are already preset.)




## qPCR Cycler Programs for Venor® GeM qOneStep Kit

- For assigning the names and/or another task (negative control, unknown, standard) to the wells switch to the tab „Assign Targets and Samples“.



- The cycling program is uploaded automatically. Click on start for PCR after loading the samples in the tray.



## qPCR Cycler Programs for Venor® GeM qOneStep Kit

---

### Programming of Rotor-Gene® 6000 (5-plex)

#### Program Step 1: Pre-incubation

Setting	Hold
Hold Temperature	95 °C
Hold Time	3 min 0 sec

**Please check the correct settings for the filter combination:**  
**green filter (510): Mollicutes**  
**yellow filter (555): Internal Control**

#### Program Step 2: Amplification

Setting	Cycling
Cycles	45
Denaturation	95°C for 5 sec
Annealing	55°C for 30 sec → <b>acquiring to Cycling A (green and yellow)</b>
Elongation	60°C for 45 sec
Gain setting	automatic (auto Gain)
Slope Correct	activated
Ignore First	deactivated

#### Result Reading:

- Open the menu *Analysis*
- Select *Quantitation*
- Check the required filter set (green and yellow) according to the following table and start data analysis by double click.
- The following windows will appear:  
*Quantitation Analysis - Cycling A (green or yellow)*  
*Quant. Results - Cycling A (green or yellow)*  
*Standard Curve - Cycling A (green or yellow)*
- In window *Quantitation Analysis*, select first *linear scale* and then *slope correct*  
Threshold setup (not applicable if a standard curve was included in parallel and auto threshold was selected)
  - In window *CT Calculation* set the threshold value to 0-1
  - Pull the threshold line into the graph. Adapt the threshold line to the initial linear section of the positive control reaction.
- The Ct-values can be taken from the window *Quant. Results*.

**qPCR Cycler Programs for Venor® GeM qOneStep Kit**

---

**Programming the Mx3005P®**

- Go to the setup menu, click on “Plate Setup”, check all positions which apply
- Click on “Collect Fluorescence Data” and check FAM™ and HEX™
- Corresponding to the basic settings the “Reference Dye” function should be deactivated
- Specify the type of sample (negative or positive control, sample, standard) at “well type”
- Edit the temperature profile at “Thermal Profile Design”:  
Segment 1: 95 °C for 3 min  
Segment 2:  
Denaturing 95 °C for 30 sec  
Annealing 55 °C for 30 sec & data collection end  
Extension 60 °C for 45 sec  
45 cycles
- at menu “Run Status” select “Run” and start the cycler by pushing “Start”

**Analysis of raw data:**

- In the window “Analysis” tab on “Analysis Selection / Setup” to analyze the marked positions
- Ensure that in window „algorithm enhancement“ all options are activated:  
Amplification-based threshold  
Adaptive baseline  
Moving average
- Click on “Results” and “Amplification Plots” for an automatic threshold
- Read the Ct values at “Text Report”

## qPCR Cycler Programs for Venor® GeM qOneStep Kit

### Programming the LightCycler® 480

#### Program 1: Pre-incubation

Cycles	1
Analysis Mode	None
<b>Temperature Targets</b>	<b>Segment 1</b>
Target Temperature [°C]	95
Incubation time [min]	3:00
Temperature Transition Rate [°C/s]	4.4
Secondary Target Temperature [°C]	0
Step Size [°C]	0.0
Step Delay [Cycles]	0
Acquisition Mode	None

#### Program 2: Amplification

Cycles	45		
Analysis Mode	Quantification		
<b>Temperature Targets</b>	<b>Segment 1</b>	<b>Segment 2</b>	<b>Segment 3</b>
Target Temperature [°C]	95	55	60
Incubation time [s]	30	30	45
Temperature Transition Rate [°C/s]	4.4	2.2	4.4
Secondary Target Temperature [°C]	0	0	0
Step Size [°C]	0.0	0.0	0.0
Step Delay [Cycles]	0	0	0
Acquisition Mode	None	Single	None

#### Program 3: Cooling

Cycles	1
Analysis Mode	None
<b>Temperature Targets</b>	<b>Segment 1</b>
Target Temperature [°C]	40
Incubation time [s]	30
Temperature Transition Rate [°C/s]	2.2
Secondary Target Temperature [°C]	0
Step Size [°C]	0.0
Step Delay [Cycles]	0
Acquisition Mode	None

**Before starting the LC480, make sure that the filter setting is correct:**

LightCycler® 480	Mollicutes	Internal Control
Instrument I	533 nm	568 nm
Instrument II	510 nm	580 nm




## qPCR Cycler Programs for Venor® GeM qOneStep Kit

### Programming the Biorad CFX 96™

#### Run Setup Protocol Tab


Click **Create New** to open the *Protocol Editor* to create a new protocol.

Select any step in either the graphical or text display — the step becomes highlighted in blue. Click the temperature or well time to directly edit the value.



#### Run Setup Plate Tab

Click **Create New** to open the *Plate Editor* to create a new plate.



Use the **Scan Mode** dropdown menu in the Plate Editor toolbar to designate the data acquisition mode to be used during the run.

Important!!! Select the *All Channels* mode.

## qPCR Cycler Programs for Venor<sup>®</sup>GeM qOneStep Kit

---

Click **Select Fluorophores** to indicate the fluorophores that will be used in the run. Choose FAM<sup>™</sup> for the detection of mycoplasma amplification and HEX<sup>™</sup> for monitoring the amplification of the internal control. Within the plate diagram, select the wells to load.

### Quantification Tab

The amplification chart displays traces of the relative fluorescence collected from each well at every cycle of the run.

Choose the fluorophore data you want to display by clicking the fluorophore checkboxes located under the amplification chart. Select FAM<sup>™</sup> to display data of mycoplasma detection and select HEX<sup>™</sup> to display internal control amplification data.

### Data Analysis

The software uses two modes for quantification cycle determination. Select Settings from the menu bar and select *Baseline Subtracted Curve Fit* as baseline setting and *Single Threshold* mode as Cq Determination Mode. In the Single Threshold mode, click and drag the threshold line to manually position the line. Adapt the threshold line to the initial linear section of the positive control reaction.

Samples showing no C<sub>t</sub>-value can be considered as negative.

### Trademarks

LightCycler is a registered trademark of a member of the Roche Group. ABI Prism is a registered trademark of Applied Biosystems Corporation or its subsidiaries in the US and certain other countries. CFX96 Touch is a trademark of Bio-Rad Laboratories, Inc. Rotor-Gene is a registered trademark of the Qiagen Group. Mx3005P is a registered trademark of Stratagene, an Agilent Technologies company. StepOne and StepOnePlus are trademarks of Applied Biosystems or its subsidiaries in the U.S. and/or certain other countries. FAM and HEX are trademarks of Applied Biosystems Corp. or its subsidiaries.

#### Minerva Biolabs GmbH

Schkopauer Ring 13 · D-12681 Berlin  
Tel. +49 (0)30 2000 437-0 · Fax +49 (0)30 2000 437-9  
info@minerva-biolabs.com · www.minerva-biolabs.com

#### Minerva Biolabs Inc.

1 Jill Ct., Building 16, Unit 10 · Hillsborough, NJ 08844 · USA  
Phone 1-908-524-4661  
info@minervabiolabs.us · www.minervabiolabs.us

© 2019 Minerva Biolabs GmbH  
Number: TN24.03EN  
Date of Release 20.06.2019  
Page 10 / 11

**qPCR Cycler Programs for Venor® GeM qOneStep Kit**

---

**Minerva Biolabs GmbH**

Schkopauer Ring 13 . D-12681 Berlin  
Tel. +49 (0)30 2000 437-0 . Fax +49 (0)30 2000 437-9  
info@minerva-biolabs.com . www.minerva-biolabs.com

**Minerva Biolabs Inc.**

1 Jill Ct., Building 16, Unit 10 . Hillsborough, NJ 08844 . USA  
Phone 1-908-524-4661  
info@minervabiolabs.us . www.minervabiolabs.us

© 2019 Minerva Biolabs GmbH  
Number: TN24.03EN  
Date of Release 20.06.2019  
Page 11 / 11