

Product Information



Internal Control extra

Supplement for Venor®GeM Classic and Venor®GeM qEP Kits

Component	Quantity	Cap color
	Cat. No. 11-9905/ 11-1905	
Internal Control extra	4 vials, lyophilized	yellow
PCR Grade Water	1 vial, 2.0 ml	white

Instructions for Use

The Internal Control extra is used exclusively in combination with Venor®GeM Classic or Venor®GeM qEP kits for the detection of mycoplasma by conventional PCR or qPCR, respectively. This product is required to carry out process validation and EP 2.6.7 / JP G3 compliant testing. Please refer to the Instructions for Use of Venor®GeM Classic or Venor®GeM qEP Kit for information about the correct rehydration and application of the product.

Storage

The lyophilized product must be stored at +2 to +8 °C. The rehydrated Internal Control extra must be stored at ≤ -18 °C. Avoid repeated freezing and thawing.

Related Products

Sample Preparation

56-1010/-1050/-1200	Venor®GeM Sample Preparation Kit	10/50/200 extractions
---------------------	----------------------------------	-----------------------

Contamination Control Kits for conventional PCR

11-1025/-1050/-1100/-1250	Venor®GeM Classic Mycoplasma Detection Kit	25/50/100/250 reactions
11-7024/-7048/-7096/-7240	Venor®GeM Advance Mycoplasma Detection Kit	24/48/96/240 reactions
11-8025/-8050/-8100/-8250	Venor®GeM OneStep Mycoplasma Detection Kit	25/50/100/250 reactions
12-1025/-1050/-1100/-1250	Onar® Bacteria Detection Kit	25/50/100/250 reactions

Contamination Control Kits for qPCR

11-9025/-9100/-9250	Venor®GeM qEP Mycoplasma Detection Kit	25/100/250 reactions
11-91025/-91100/-91250	Venor®GeM qOneStep Mycoplasma Detection Kit	25/100/250 reactions

Please visit our website for more information about our products.

Document Version: 01 | Document Release 30.09.19

Germany:

Minerva Biolabs GmbH

Schkopauer Ring 13 | 12681 Berlin | Germany

Tel.: +49-30-2000437-0 | Fax: +49-30-2000437-9

E-mail: order@minerva-biolabs.com

Internet: www.minerva-biolabs.com

USA / Canada:

Minerva Biolabs Inc.

1 Jill Ct., Building 16, Unit 10

Hillsborough | NJ 08844 | USA

Phone: 1-908-524-4661

order@minervabiolabs.us | www.minervabiolabs.us